
WM 66610a 11/2017 DE, EN, FR

CE 0197

Löwenstein Medical
Technology GmbH + Co. KG
Kronsaalweg 40
22525 Hamburg, Germany
T: +49 40 54702-0
F: +49 40 54702-461
www.loewensteinmedical.de

OXYnator supra

Nasal-/Oralbrille

Nasal/oral cannula

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WM 66610a


1


2

①

②

1 Bedienung

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1 Operation

Look at the diagrams for how to put on the nasal/oral cannula:

- 1 Arrange the nasal/oral cannula as shown
- 2 Put the nasal/oral cannula in position

Proceed in reverse sequence to remove the nasal/oral cannula.

2 Intended use

The nasal/oral cannula is intended for the nasal/oral supply of oxygen to individuals weighing over 30 kg. The nasal/oral cannula is intended for use on an individual patient and not for a change of patient. When used in clinical facilities without a change of patient, the nasal/oral cannula is replaced after 72 hours. In the domestic environment, the nasal/oral cannula is intended for a single patient and a maximum service life of 30 days.

3 Description of function

When oxygen is supplied and the nasal/oral cannula is attached in the area of the top lip, it generates a cloud of oxygen which can be inhaled through the nose or mouth. This method is scientifically proven; there is evidence of an adequate therapeutic effect even if the nasal passages are blocked.

4 Safety

⚠ WARNING

Risk of fire from the use of oxygen in combination with flammable substances
Supplying oxygen without a special safety device can lead to fire and injure people.
⇒ Do not operate the nasal/oral cannula in the vicinity of excessive heat, naked flame or objects at risk of fire.

⚠ WARNING

Risk of infection if the nasal/oral cannula is reused!
If the nasal/oral cannula is used for too long, infections may be transmitted.
⇒ Regularly replace the used nasal/oral cannula with a new one before unhygienic contamination occurs.

⚠ WARNING

Risk of injury if the nasal/oral cannula malfunctions
A damaged nasal/oral cannula or damaged individual parts may injure the patient.
⇒ Only use undamaged nasal/oral cannulas.
⇒ Do not crush or kink tubes, cannulas, and connections.
⇒ Keep away objects on which the tubes may get caught.

Read and follow the complete Instructions for Use before using the nasal/oral cannula.

5 Product description

5.1 Overview

A diagram of the individual parts can be found on the title page.

- 1. Connection to oxygen supply
- 2. Attachment sleeve
- 3. Double patient circuit
- 4. Base plate with ball joint tube connections

6 Cleaning

It is not possible to clean or disinfect the entire nasal/oral cannula reliably.

The manufacturer recommends cleaning just the outside of the nasal/oral cannula without its base plate under running water or by wiping over with a solution of alcohol-based disinfectant.

You can clean the base plate without the ball joint tube connections as outlined below.

1. Pull both ball joint tube connections off the base plate.
2. Wash the base plate in warm water and a mild detergent.
3. Disinfect the base plate in a solution of alcohol-based disinfectant.
4. Rinse the base plate with clean water.
5. Allow the base plate to air-dry.
6. Subject the base plate to a visual inspection.
7. Push both ball joint tube connections onto the base plate.

7 Disposal

In a domestic environment, you can dispose of all parts in domestic waste.

In a clinical environment, dispose of all parts in line with clinic regulations.

8 Technical specifications

Product class to Directive 93/42/EEC	Ila
Flow rate	0 l/min to 3 l/min
Storage temperature	-25 °C to +70 °C
Operating temperature	+5 °C to +30 °C

9 Materials

Part of the nasal/oral cannula	Material
Patient circuit and connections	PVC (polyvinyl chloride)
Attachment sleeve, ball joints and base plate	PP (polypropylene)

10 Declaration of conformity

The manufacturer Löwenstein Medical Technology GmbH + Co. KG (Kronsaalsweg 40, 22525 Hamburg, Germany) hereby declares that the product complies with the relevant regulations of Directive 93/42/EEC governing medical devices. The unabridged text of the Declaration of Conformity can be found on the manufacturer’s website.

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