

FreeCheckTest OxyCare GmbH

PatientData

Name : Mustermann1
ID : 1
Physician :
Note1 :

Gender : male
Height : 0 in

Age : 0
Weight : 0 lb

Note2 :

Start Date: 10/09/08

Start Time : 19:22:57

End Date: 10/09/08

End Time : 07:58:29

Comments:



Event Data	Spo2	Pulse			
Total Events.	4	6	Spo2 > 94	36875	times
Time In Events.(Min)	1	175	Spo2 < 88	41	times
Avg.Event Dur.(Sec)	11165	7443	Pulse > 150	0	times
Index.(1/hr)	0.3	0.5	Pulse < 60	1357	times
Avg.Event Max.	95	112	Diagnose:		
Avg.Event Min.	85	56			
Max.Event offset.	15	39			
Max Time.Event.(Min)	0	106			
Max Dur.	99	139			
Min Dur.	78	36			
Avg Dur.	95	66			

Desaturation Event: drop in Spo2 by at least 4% for a minimum duration of 10 seconds.



Pulse Event Change in rate by at least 6 bpm for a minimum duration of 8 seconds.

Graphic Summary

Spo2 (10% per division)



Pulse Rate (20BPM per division)



20:37:10

23:07:14

01:41:32

04:08:54

06:37:05