



# SISS BABYCONTROL<sup>®</sup> M

 **MASIMO SET<sup>®</sup>**

**OXYCARE<sup>®</sup>**  
**Medical :: Group GmbH**  
FON +49(0)421 48 99 66 • [www.oxygene.eu](http://www.oxygene.eu)

MEMBER  
OF



**SCHULTE-ELEKTRONIK GMBH**  
MEDICAL TECHNOLOGY

# SISS BABYCONTROL® M



RESPIRATION  HEART RATE  OXYGEN SATURATION SpO<sub>2</sub>

**WILDDESIGN®** and 

**SISS BABYCONTROL® M** is the latest handheld patient monitor, developed in collaboration with WILDESIGN, a specialist in medical design, and MASIMO, the well known manufacturer of SpO<sub>2</sub> technologies.

## Mobility

With a weight of only 450 grams the monitor operates with either regular or rechargeable batteries. Optimal mobility during monitoring is given.

## Clear and easy

The graphical LCD-Display and the clearly arranged control elements allow an intuitive operation. Colored LEDs clearly show the status of the monitor. The **SISS BABYCONTROL® M** gives optical and acoustical alarm depending on the event.

## Easily removable MicroSD memory card

Despite its small size the **SISS BABYCONTROL® M** measures and stores complex procedures automatically. The MicroSD memory card can be easily removed from the device. To evaluate the stored data with the help of a functional software it can be uploaded on a PC.

## Technical Data

Display	Vital signs via LCD / Alarms via LED
SpO <sub>2</sub> Connector	MASIMO SET® technology
Respiration Connector	Special Respiration Sensor
ECG Connector	2/3 Disposable ECG Electrodes
Memory	2 GB MicroSD memory card
Weight / Dimensions	450 grams / 150 x 80 x 55 mm
Power supply	4 x AA Battery / Akku



**SCHULTE-ELEKTRONIK GMBH**  
MEDICAL TECHNOLOGY

Schörenbergstr. 20 • 59939 Olsberg • Germany  
Tel. +49 (0)2962 9707 0 • Fax +49 (0)2962 9707 15  
www.schulte-elektronik.de • eMail: info@schulte-elektronik.de