

## Features and Benefits

### F&P Pilairo Q Adjustable Headgear – 400HC326

Adjustable headgear fits a greater range of head sizes and enables sizing adjustment to patient preference.

### F&P AirPillow™ Seal – 400HC125


Inflatable AirPillow Seal self-inflates in and around the nose. Can move and flex without compromising the seal.

### F&P Pilairo Q Nasal Pillows Mask (no headgear) – 400HC565

Low profile and lightweight Easy Frame for a clear line of sight



Embedded VELCRO®\*\* Tabs enable fast fitting and removal.




Ball-and-Socket Elbow rotates for freedom of movement while reducing drag from the CPAP tube.



### F&P Pilairo Q Stretchwise™ Headgear – 400HC322

Minimal Stretchwise™ Headgear automatically adjusts to patients for a quick and easy fit.

### F&P Pilairo Q with Adjustable & Stretchwise Headgear – 400421



Mask stabiliser for additional stability



### F&P Pilairo Q Elbow Cover & Diffuser (10 pack) – 400HC229

Q Cover for quieter use<sup>1</sup> and minimal draft.

Lightweight flexible tube.

### Female Swivel – 400HC224

Swivel with a quiet and smooth rotation.

\*VELCRO is a registered trade mark of Velcro Industries B.V.

1. The Sound Pressure Level of the mask with the Q Cover is 23.6 dBA, with uncertainty 2.5 dBA.