

# FlexiFit™ 431

FULL FACE MASK

## Specification Sheet

A Full Face Delivery Mask




*Incorporating the FlexiFit™ Technology with auto-contouring properties, the FlexiFit™ 431 Full Face Mask comes with an Under-Chin fit and highly contoured Silicone Seal for maximum comfort and stability.*

**SleepStyle™**  
SLEEP PERFORMANCE FOR  
AN ENERGIZED LIFESTYLE


**Fisher & Paykel**  
HEALTHCARE

# Full Face Delivery - *FlexiFit™ 431*




*The FlexiFit™ 431 Full Face Mask offers a comfortable alternative to nasal or oral delivery and comes as an Out-of-Box Solution fully assembled and ready for use.*

## UNIQUE MASK FEATURES FOR GREATER FREEDOM, COMFORT AND COMPLIANCE




### FlexiFit™ Technology Auto-contouring fit

*With Air Cushion*  
FlexiFit™ Technology automatically adjusts to a wide range of nasal and facial contours with simple one-step ease of fitting and no complicated adjustments.



### Glider™ Strap Anti-leak system

The sliding Glider™ Strap provides freedom of movement while maintaining a seal. Mask force remains equal when the head rotates from side to side, minimizing the occurrence of leaks.




### Under-Chin Design Extra stability

An Under-Chin fit provides stability by providing a constant point of reference and maintains the seal even with movement of the jaw during the night.



### Efficient Bias Diffuser Quiet usage

*Maintenance Free*  
This advanced air diffuser system prevents air from blowing onto the sleeping partner and decreases noise that could potentially disrupt sleep.




### Premium Frosted Silicone Seal Optimizing seal technology

This seal technology has a soft, smooth feel to slide onto the face, and its frosted texture optimizes the mask seal.



### Crown Strap Extra stability

Complements the Under-Chin Design to maximize stability and eliminate the need for a chin strap.



### Easy-Release Strap One-step removal

Simply pull the Easy-Release Cord for one-step removal of headgear. No ongoing headgear adjustments required.

## CONTENTS

- Complete mask with headgear assembled and ready for use
- Alternate size Silicone Seals included (S,M,L)
- User friendly Fitting Guide

REORDER CODE: HC431

FOR MORE INFORMATION, PLEASE CONTACT YOUR LOCAL FISHER AND PAYKEL HEALTHCARE REPRESENTATIVE:

**International** PO Box 14 348, Panmure, Auckland 1134, New Zealand Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: [info@fphcare.com](mailto:info@fphcare.com) Web: [www.fphcare.com](http://www.fphcare.com)

**USA** Tel: 1800 446 3908 +1 949 470 3900 Fax: +1 949 470 3933 **Australia** Tel: +61 3 9879 5022 Fax: +61 3 9879 5232 **Austria** Tel: 0800 29 31 23 Fax: 0800 29 31 22 **Benelux** Tel: +31 78 644 0924

Fax: +31 78 644 0914 **China** Tel: +86 20 8363 4286 Fax: +86 20 8363 3595 **France** Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 **Germany** Tel: +49 7182 93777 0 Fax: +49 7182 93777 99

**India** Tel: +91 80 5123 6041 Fax: +91 80 5123 6044 **Irish Republic** Tel: 1800 409 011 **Italy** Tel: +39 06 7839 2939 Fax: +39 06 7814 7709 **Japan** Tel: +81 3 3661 7205 Fax: +81 3 3661 7206

**Spain** Tel: +34 902 013 346 Fax: +34 902 013 379 **Switzerland** Tel: 0800 83 47 63 Fax: 0800 83 47 54 **UK (EU Authorised Representative)** Tel: +44 1628 626 136 Fax: +44 1628 626 146

**Nthn Ireland** Tel: 0800 132 189

