

# PHILIPS

Over-the-nose  
nasal mask

Pico

Traditional nasal mask

Tube-in-front design

Adjust with headgear clips

HH1032/01



## Value-added performance

A truly appealing traditional nasal mask that makes minimal contact with the face and offers simplicity, comfort and seal.

### Compact design

- Pico comfortably fits more than 98% of users\*
- The lightest and smallest traditional nasal mask \*

### Performance and comfort

- Provides for quick and easy fitting and adjustment

## Highlights

### Compact design for a natural fit




The 'tip-of-the-nose' cushion is designed to sit below your nose bridge, making minimal contact with your face while creating a superior seal. \* Data on file, Internal Marketing Survey 2013

### Minimal design, superior seal



The design and headgear clips make it easier to use and clean.

### The lightest, smallest choice



Pico makes minimal contact with the face, providing the freedom to watch television or read while wearing it. \* Data on file, Internal Marketing Survey 2013

## Specifications

### Service

Warranty: 90-days

### Specifications

Operating pressure: 4 -20 cm H2O

Sound pressure level: 19 dBA

N° of parts: Five

Cushion material: Silicone/polycarbonate

Headgear material: UBL/Urethane

foam/nylon/Lycra

**Cleaning instructions: mask:** Hand wash with mild detergent, Rinse thoroughly, air dry

**Cleaning: headgear, tubing:** Hand wash with mild detergent, Rinse thoroughly, hang to dry

**Frame/Elbow/Swivel material:** Polycarbonate

**Talon clips material:** Thermoplastic polyester elastomer

### Intent of use

To be used on patients <gt/>30kg: with CPAP prescription

### Light and small

**Easy fitting and adjustment:** Quiet exhalation

