

PHILIPS

PEP

Threshold

Positive Expiratory Pressure
device

For airway clearance

Helps remove lung secretions

HH1333



Pressure that's good for you

With a flow-dependent, one-way valve, threshold devices maintain a constant pressure, regardless of the patients airflow.

Convenient design

- Reliable and easy to clean
- Works in any position

Mouthpiece/mask option

- Mouthpiece/mask option

Intended use

- Intended use of PEP

Highlights

Mouthpiece/mask option

Mouthpiece/mask option for flexibility

Reliable and easy to clean

Threshold is made of durable, high-impact acrylic.

Works in any position

Effective therapy in any position

Intended use

Threshold PEP expiratory pressure trainer for can be used for airway clearance, bronchiole hygiene and chest physical therapy.

Specifications

Product details

Warranty: Replace after 1 year

Intended Use: As directed by your physician

Cleaning: Wash in warm soapy water after each use. Flush all parts with clean water.

Shake out excess water and air dry. DO NOT BOIL OR HEAT.

Includes: training diary

