



JOYCEeasy

... ingenious simplicity.

Nasal or Full Face. Simply easy.

It could not be any simpler. Everything about our ventilation mask JOYCEeasy is designed for easy use – from simple mask adjustments directly on the patient's face to intuitive handling. Because the mask fits so well, JOYCEeasy also encourages better therapy compliance by minimizing leaks and pressure points. The quiet, diffuse exhalation system ensures a peaceful bedroom atmosphere.

Ingenious simplicity everywhere:

- Unlimited freedom of movement provided by the ball-and-socket joint **1**
- Quiet and diffuse exhalation system **2**
- Mask can be put on and taken off with just one hand
- Color-coded hook-and-loop tabs **3**
- Extremely robust plastic
- Quick-release cord (only Full Face)

JOYCEeasy next FF




JOYCEeasy Full Face



JOYCEeasy nasal

**NEW
DESIGN**



JOYCEeasy next FF

JOYCEeasy Full Face

... the perfect duo.

The right fit for everyone.

We heard you.

And we reworked JOYCEeasy Full Face to satisfy your practical needs.


JOYCEeasy next FF now features a completely redesigned large mask cushion that has a comfortable and firm fit even at high pressures.

The headgear too remains securely in place at high pressures. Thanks to the right material mix and ergonomic design, the headgear clips can be easily loosened from the mask.

Your benefit – The right mask for everyone:

JOYCEeasy next FF = normal to obese

JOYCEeasy Full Face = slender to gaunt



JOYCEeasy nasal

Thanks to your feedback, we have made more improvements to the JOYCEeasy mask cushion.

The wider and softer cushion sits lower across the bridge of the nose and easily accommodates larger noses comfortably. The flatter and more flexible upper lip area lies more gently on the face.

Now even more patients can benefit from the enlarged mask cushion.

Don't worry! The size of the cushion has not been changed in Size S, which still fits small faces.

Accessories

Ball-and-socket set, vented nasal retaining ring, elbow, rotating sleeve
WM 15798



Ball-and-socket set, vented Full Face (JOYCEeasy next FF) retaining ring, elbow, anti-asphyxia valve, rotating sleeve
WM 15944



Ball-and-socket set, vented Full Face (JOYCEeasy Full Face) retaining ring, elbow, anti-asphyxia valve, rotating sleeve
WM 15897



Quick-release cord
WM 26890 (Full Face)



Endoscope adapter set JOYCEeasy next FF NV / JOYCEeasy Full Face NV retaining ring NV, endoscope elbow NV
WM 15864




Elbow set JOYCEeasy next FF NV / JOYCEeasy Full Face NV retaining ring NV, elbow NV
WM 25269




Replacement parts


WEINMANNstrap for JOYCEeasy nasal including headgear clips
WM 25505




WEINMANNstrap for JOYCEeasy Full Face including headgear clips
WM 25606




WEINMANNstrap for JOYCEeasy next FF including headgear clips
WM 25105



Mask cushion silicone
Size S WM 26803 (Nasal)
WM 26853 (Full Face)
WM 25104 (next FF)
Size M WM 26813 (Nasal)
WM 26863 (Full Face)
WM 25114 (next FF)
Size L WM 26823 (Nasal)
WM 26873 (Full Face)
WM 25124 (next FF)



Forehead cushion silicone
WM 26801 (Nasal and Full Face)



Product name	Size S	Size M	Size L
JOYCEeasy Nasal	WM 25510	WM 25520	WM 25530
JOYCEeasy Full Face	WM 25910	WM 25920	WM 25930
JOYCEeasy next FF	WM 25109	WM 25119	WM 25129
JOYCEeasy next FF NV	WM 25159	WM 25169	WM 25179

For more information on our therapy solutions, accessories and mask systems please visit loewensteinmedical.de

CE 0197



p10105en1705

© Copyright protected. Copying or reproduction of any sort requires the explicit approval of Löwenstein Medical Technology.

Löwenstein Medical UK
1 E-Centre, Easthampstead Road
RG12 1 NF Bracknell
United Kingdom
T: +44 1344 830023
info@loewensteinmedical.uk
www.loewensteinmedical.uk

Löwenstein Medical Technology
Kronsaalsweg 40, 22525 Hamburg
Germany
T: +49 40 54702-0
F: +49 40 54702-461
info@loewensteinmedical.de
www.loewensteinmedical.de



JOYCEeasy

Easy to Fit, Easy to Use.

